

Migratek Legacy PACS Connector Status 8/22/2011

	dicomBLAST	Semperdata		
Legacy System	Status	Status	Availability	Notes
Generic DICOM (CMOVE)	Done	Done	Now	
Generic DICOM (File Tree)	Done	Done	Now	Configurable
Agfa Heartlab Cardiology	Done	Done	Now	
Agfa IMPAX	Done	Done	Now	
Amicas PACS	Not planned	Done	Now	
Cerner ProVision	Not planned	Done	Now	
DR Systems	Not planned	Done	Now	
eRad PACS	Not planned	Done	Now	
Fuji Synapse	Done	Done	Now	
GE Centricity 2.x	Done	Done	Now	
GE Centricity 3.x	Not planned	Done	Now	
GEMNET Cardiology	Done	Done	Now	
Kodak Archive Manager (Cemax-Icon)	Done	Not planned	On request	
McKesson HMI v4.x and earlier/ Camtronics Vericis	Not planned	Done	Now	
Merge eMed	Done	Done	Now	
OpenMed PACS	Not planned	Done	Now	
Philips (Sectra) EasyAccess	Not planned	Done	Now	
Philips (Stentor) iSite	Not planned	Test	2011 Q3	
Philips Inturis Cardiology	Done	Done	Now	
Siemens KinetDx	Done	Done	Now	
Siemens Sienet (MagicStore)	Done	Done	Now	
Stair Systems	Done	Not planned	On request	

Note: Laitek's batch-oriented "dicomBLAST" platform is not used for new designs, but may be used for smaller projects involving older or rare legacy systems.